


ImmuSpeed™: Novel Automated Platform for Fast ELISA in Microchips

Joël Rossier, Christine Vollet, Patrick Morier
Valdone Meylan and Frédéric Reymond

DiagnoSwiss S.A., Route de l'Île au Bois 2, CH-1870 Monthey (Switzerland)
info@diagnoswiss.com


ImmuSpeed™ is an automated microfluidics-based platform, designed for running **bead-based ELISA** tests with standard immunology reagents. The instrument works with a microfluidics chip, harboring biosensor channels, in which all assay steps, from sample introduction through signal detection, are performed automatically. Programming of sequential assays in **8 parallel channels** and regeneration loops allows running **96 tests in 1-1.5 h, fully unattended**. Various assays are available by DiagnoSwiss including: **TSH, T3, T4, LH, FSH, hCG, AFP, Testosterone, Progesterone, Prolactin**.

Principle of magnetic bead based ELISA



- Samples distributed in 96-well plate and automatically processed in batch
- Possibility to run multi-menu assays according to choice of antibody

ImmuSpeed Chip: 8-microchannel biosensor



- Beads trapped within ~300 nL channel, with large pre-concentration
- >10-fold reduction in consumption of ELISA kit reagents

ImmuSpeed Instrument



Picture of ImmuSpeed™ instrument in its 4 micro-liter plate configuration


Microfluidic platform for automated ELISA comprising:

- 2, 3 or 4 positions for 96-well plates
- Slide-in ImmuSpeed-Chip holder
- Shakers for sample & mixing plates
- Robotic station for sample & bead mixing, reagent dispense and assay performance
- Multichannel amperometric reader
- Temperature control
- Integrated software ImmuSpeed-Soft with dedicated assay protocols

Footprint: 80 x 40 x 45 cm³
Operating temperature: 37°C
Feedback control of automation

Software: ImmuSpeed-Soft, integrated in the instrument, to control the robotics and to run dedicated protocols experiments and to treat the measurement data.

Process Workflow in ImmuSpeed



Top view of ImmuSpeed™ with typical plate and chip positioning

Typical procedure for sandwich-type immunoassays:

- 1 Samples of Row 1 dispensed in mixing plate
- 2 Beads of Row 1 mixed with samples in mixing plate
- 3 Conjugate mixed with samples & beads in mixing plate
- 4 Complex dispensed in microchip
- 5 Washing and substrate addition
- 6 Electrochemical detection in chip

Example of Application: TSH tests

The convenience and benefits of the ImmuSpeed platform are demonstrated here by the detection of Thyroid Stimulating Hormone (TSH) in plasma samples.


The samples and antibody-coated magnetic beads are dispensed in 96-well plates and placed in the instrument, along with a clean mixing plate and a deep plate containing the washing buffer and the enzymatic substrate solution. ImmuSpeed™ then drives the assay automatically according to the above workflow (with 25 minutes of bead-conjugate-sample incubation for each assay).

After the incubation, 20 µL of the mix are dispensed in the chip for bead capture, washing and electrochemical detection. After each series of 8 tests, the beads are washed out of the biosensor channels and the chip is regenerated using substrate solution.

Assay Characteristics:

- Throughput of 1 plate / hour, incl. chip regeneration
- LOD in the low pg/mL range
- CV < 10; ,Dynamic range: >3 orders of magnitude
- Hands-on time: 10 min (sample & bead dispense only)
- Sample volume: from 10 µL
- Calibration: standards provided for each type of assay

- Reagent shelf-life: 12 months
- Automatic, all-in-one solution (incl. mixing, incubation, washing, detection, and ev. dilution)



Calibration curve (A) and standard deviation (B) for the detection of TSH in plasma samples using antibody-coated superparamagnetic beads (MyOne Dynabeads, Invitrogen, at 0.1 mg/mL), ALP conjugated detection antibodies (bioMérieux at 1:100), washing buffer and substrate solution (p-aminophenyl phosphate, DiagnoSwiss, at 5 mM pH 9.8)
(Insert: color plate representation for repeated assays as dispensed in the microtiter plate)

Conclusion: ImmuSpeed: an automated platform for rapid ELISA in small volume

ImmuSpeed™ is a robotized microfluidic analytical platform dedicated to fast ELISA in low volumes. By dramatically reducing hands-on time and analysis, significantly decreasing sample/reagent consumption and cost, this novel biosensor system enables to perform multi-menu immunoassays with easy-to-use robotics and laboratory infrastructure.

References

- J. S. Rossier et al., J. Assoc. Lab. Automation, 2008, 13(6), 322-329
- J. S. Rossier and H. H. Girault, Lab Chip, 2001, 1, 153
- D. Hoegger et al., Anal. Bioanal. Chem., 2007, 387, 267
- P. Morier et al., Electrophoresis, 2004, 25, 3761

Acknowledgements

DiagnoSwiss thanks the European Commission for financial support of part of this work (Grants 037953 "NeuroTAS", 201333 "DECanBio" and 227243 "Nanobe").

